

Immedia SmartMove

SUPPORT PAD WITH HANDLES

Immedia SmartMove acts as an extension of the arms, allowing the caregiver to offer support and assistance whilst maintaining a good ergonomic position and grip.

SmartMove pad can be used to support and help transfer from sitting to standing, standing to sitting, to move higher up in bed or move back/forward in a

chair. SmartMove provides the caregiver with a place to grip and to maintain a good working position. The inside has a non-slip material to help prevent sliding out of position during the transfer, whilst the outside has a low friction side to aid application and removal.



Immedia SmartMove has a non-slip material to help prevent sliding out of position.

INTENDED USE

- Support stand up / sit down
- Support moving higher up in bed
- Moving back / forward in chair / wheelchair
- Support transfer between bed, wheelchair, chair, etc

| DESCRIPTION | WIDTH X LENGTH | | ITEM NO. | SPECIFICATIONS |
|---|----------------------|---------------------------|----------|---|
| IMMEDIA SMARTMOVE ~ 4Handles | CM 21 x 60 | INCH 8.3 x 23.6 | IM426 | OUTSIDE: Nylon with polyurethane coating INSIDE: Polypropylene HANDLES: Polyester |
| WASHING INSTRUCTIONS Do not use fabric softener or bleach. <div> </div> | | | | |
| | | | | |